

DELTAICE_{2L}

Smart Passive Blood Storage Unit

REMOTE ALERTS

Bluetooth and Cellular alerts warn you and your command of temperature excursions

TO BE FDA LISTED

Demonstrate safety simply by showing our FDA listing to your blood supplier

EXTREME HEAT ENVIRONMENTS

Stays cold 12 hours even in uncooled vehicles

ALL ORIENTATIONS

Temperature maintained in any orientation

RUGGED DESIGN

Insulation protected from drops and punctures

KNOW YOUR STORAGE TIME

Screen and app informs remaining time before thermal storage melts

RECHARGEABLE BATTERY


Recharge with USB type C

LONG STORAGE TIME

Stays cold 30 hours in typical cooled vehicles and records temperature for 72 hours before recharging



911 Tactical Medicine
www.911tacmed.com
270-871-2953



GENERAL**Size**

11.5"L x 9.7"W x 9.5"H

Weight

10.5 lbs unloaded

FDA Classification

Class II. Code KSE

ENVIRONMENTAL**Temperature Ranges:****Operating**

2° to 60°C (36°F to 140°F)

Storage

-20° to 70°C (-4°F to 158°F)

Ingress Protection

IP65, washdown capable

STORAGE**Duration**

48 hours in 25°C (77°F)

24 hours in 44°C (111.2°F)

12 hours in 60°C (140°F)

Internal Volume

2 L

Unit capacity

Three 450 ml bags

PREPARATION**Freezer Space**

4.7 L

Freeze Time

8 hrs at -9.5°C (15°F)

Prep Time15 minutes in 25°C (77°F)
ambient**TEMPERATURE RECORDING****Frequency**

15 minutes

Duration

72 days

Accuracy

± 1°C

ReportingCellular and Bluetooth
to Application